## **DIALUNOX**

# Use of a microfluidic disc in Point-of-Care applications

Workshop "Medical Microtechnology" February 22, 2023 Benjamin Albrecht – Dialunox GmbH

Use of a microfluidic disc in Point-of-Care applications

#### **Content**

- · Company profile
- Development of the discreader platform
- From the application to product requirements
- · Different fields of application

## DIALUNOX @ a glance



- Located in the south of Germany
- Founded 1997
- > 80 employees
- Since 2010 part of QIAGEN

- CDMO (Contract Development and Manufacturing Organisation)
- Engineering Services for optical measurement instruments
- Contract manufacturer for instruments







- One-Stop-Shop for customers from all industries
- Broad expertise in engineering (Software, mechanics, optics, electronics)
- State of the art production facility
- ISO 9001 & EN13485 certified

DIALUNOX

Workshop "Medical Microtechnology", B. Albrecht

February 22, 2023 3

Use of a microfluidic disc in Point-of-Care applications

## **DLX Product Portfolio**

#### **OEM Instrument Supply**



Sold > 20.000 units

- 5 different instrument platforms

#### **Optical Detector**



Built > 40.000 Units

- Unique highly sensitive fluorescence detector

#### **Engineering Service**



#### 25yrs experience

 Starting point of all business/ customization tasks Use of a microfluidic disc in Point-of-Care applications

## Development of the discreader platform

- BMBF and EU funded project, started in 2012
- Record player for a fast diagnosis
- For Immunoassays as well as PCR







Workshop "Medical Microtechnology", B. Albrecht February 22, 2023

Use of a microfluidic disc in Point-of-Care applications

## Development of the discreader platform

- Development Project started in 2017
- Triangle project
- Initial plan: MRSA -> Covid
- Increase of test portfolio



#### From the Application to product requirements

Sample Prep

DNA Extraction

**PCR** 

Result

Different type of samples

- \* swab, skin, liquid
- \* Weight in
- \* Homogenize
  - -> Spinning

Extract DNA from cellular components (Lysis)

- mechanical
- enzymatic
  - -> Positioning

- \* Denaturation
- \* Annealing
- \* Elongation
- -> 30-40 cycles
- -> Temperature control
- \* Calculation of C<sub>⊤</sub> value
- \* Display of the Amplification curves
- \* Connection to LIS
- -> Fluorescence detection

**DIALUNOX** 


Workshop "Medical Microtechnology", B. Albrecht February 22,

Use of a microfluidic disc in Point-of-Care applications

## From the Application to product requirements

- Spinning -> Up to +/- 60Hz to create the correct centrifugal force
  - · Pumping liquid from one chamber to the next
  - Rehydration of reagent in the stickpack and release
- **Heating and Cooling Unit** 
  - · Different zones for independent temperature control
  - · Fast heating and cooling rates for good **TimeToResult**
  - · Good heat transfer from heating plate to liquid
- Fluorecencs detection
  - 4 channel Multiplexing
  - · Fast measurement frequency





## Different fields of application





#### Spindiag (2016)

Molecular diagnostics at point-of-care for infectious diseases





#### Dermagnostix (2020)

Identification of inflammatory skin diseases through molecular diagnostics

## Endress+Hauser 🖽



#### Endress+Hauser BioSense (2021)

Rapid on-site molecular analyses for bacterial or viral contamination in water and beverages as well as genetic changes in food

**DIALUNOX** 

Workshop "Medical Microtechnology", B. Albrecht February 22, 2023



info@dialunox.com



Benjamin Albrecht +49 7771 / 9166 - 236 benjamin.albrecht@dialunox.com